



Quick Setup Guide

CleveMed SleepView™ Home Sleep Test

Belt Placement



Nasal Cannula Placement



Finger Probe Placement



Starting the Test

Press the “ON” button. The power light should turn green. The sensor lights will go off after 90 seconds. The power light will remain on.

Ending the Test

At the end of the test, press and hold the “OFF” button to turn off the device. If the power light has already turned off, turn it on and then off again.

When the test is finished, put everything back in the bag including the tape and documentation log, unless you’ve been given different instructions. You can find additional information at BayCareSleep.org. Scroll halfway down the page and select “Patient Instructions” to see video instructions or scan this code:



To reach us after hours: (727) 734-6716, option 3

If an emergency medical situation occurs during testing, call 911.

 **BayCare**
Sleep Disorders
Centers

BayCareSleep.org